



Bariatric Battery-Powered Lift with 4-Point Cradle

- To charge, simply plug in unit Battery does not need to be removed to be charged
- Lift will not operate while plugged in Ideal for home-care use
- Jumbo actuator pump with emergency button can switch to manual mode to lower patient down safely
- Base legs adjust easily and lock securely into open position with padded shifter handle
- Extra wide coated cradle with 360° rotation with 4-point hook up (Can fit under bed)
- Weight Capacity 600 lbs

PRODUCTS

	Description
13244	Bariatric Lift Elect Homecare Style, 1/cs

SPECIFICATIONS

	Description
Boom Height	32" - 75.5"
Base Width	45.75" (Open) 25" (Closed)
Base Height	4.5"
Base Clearance	2"
Casters	3" Front 4" Rear
Weight	132 lbs.
Carton Shipping Weight	148 lbs.
Weight Limit	600 lbs.
Warranty	Lifetime on frame
Lift Frame	Lifetime Limited
Electronics	1 yr
Battery	1 yr

PARTS

	Description
13240FP	Frame Plug, 1/ea
13244CAB	Caradle Attachment Bolt, 1/ea
13244MHB	Mast Hinge Bolt, 1/ea
13244UBA	Upper Actuator Bolt, 1/ea
13053PB	Pivot Bolt, 1/ea
13244PH	Push Handle, 1/ea
13240BATT	Rechargeable Battery, 1/ea
13240BC	Battery Cell * *Need to order 2 each, 1/ea
13240C	Charger for 13240, 1/ea
13255	External Charger, 1/ea
13240A	Actuator, 1/ea
13240HC	Hand Control, 1/ea
13244SB	Spreader Bar, 1/ea

PARTS	Description
13240SBH	Spreader Bar Handle, 1/ea
13244BAR	Base Actuating Rod, 1/ea
13244CAL	Caster with Lock, 1/ea
13244ML	Lower Mast, 1/ea
13244CPL	Caster Position, 1/ea
13240CA	Caster without Lock, 1/ea
13244B	Base
13240PC	Power Cord for 13240 charger, 1/ea
13255PC	Power Cord for 13255 external charger, 1/ea
13244RPA	Receiver Plate assy., 1/ea
13240ARB	Receiver Plate Bolt, 1/pkg
13240ARB-N	Ser # 4A1206 & Later, 1/pkg
13240AP	Spreader Locking Plate, 1/ea
13244MNB	
CPCFUSE	Fuse
13244BMB	Base to Mast Bolt
13244F	Footplate
13244CA	Non-locking Caster